

Product Datasheet

Human COL9a3 ELISA Kit (Ready to Use) (orb550306)

Description	Human Collagen Type IX Alpha 3 (COL9a3) Sandwich ELISA KitReady to Use: Pre-diluted detection reagen...
Reactivity	Human
Range	0.781-50ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q14050
Sensitivity	0.29ng/mL
Entrez	1299
Expiration Date	Please enquire.