

Product Datasheet

Human C1qA ELISA Kit (Ready to Use) (orb550328)

Description	Human Complement Component 1, Q Subcomponent A (C1qA) Sandwich ELISA KitReady to Use: Pre-diluted de...
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityAssay Length: 1-3 hoursMethod: sandwichserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P02745
Sensitivity	0.059ng/mL
Entrez	712
Expiration Date	Please enquire.