

Product Datasheet

Human C3 ELISA Kit (Ready to Use) (orb550334)

Description	Human Complement Component 3 (C3) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for...
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P01024
Sensitivity	0.061ng/mL
Entrez	718
Expiration Date	Please enquire.