

Product Datasheet

Human CRNN ELISA Kit (Ready to Use) (orb550368)

Description	Human Cornulin (CRNN) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for convenient ...
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q9UBG3
Sensitivity	0.123ng/mL
Entrez	49860
Expiration Date	Please enquire.