

Product Datasheet

Human CRIM1 ELISA Kit (Ready to Use) (orb550428)

Description	Human Cysteine Rich Transmembrane BMP Regulator 1 (CRIM1) Sandwich ELISA KitReady to Use: Pre-dilute...
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q9NZV1
Sensitivity	0.133ng/mL
Entrez	51232
Expiration Date	Please enquire.