

Product Datasheet

Human CYP11A1 ELISA Kit (Ready to Use) (orb550440)

Description	Human Cytochrome P450 11A1 (CYP11A1) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents ...
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Sensitivity	0.133ng/mL
Expiration Date	Please enquire.