

Product Datasheet

Human CYP1A1 ELISA Kit (Ready to Use) (orb550442)

Description	Human Cytochrome P450 1A1 (CYP1A1) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents fo...
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P04798
Sensitivity	0.056ng/mL
Entrez	1543
Expiration Date	Please enquire.