

Product Datasheet

Human CYP1A2 ELISA Kit (Ready to Use) (orb550443)

Description	Human Cytochrome P450 1A2 (CYP1A2) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents fo...
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P05177
Sensitivity	0.051ng/mL
Entrez	1544
Expiration Date	Please enquire.