

Product Datasheet

Human CYP1B1 ELISA Kit (Ready to Use) (orb550444)

Description	Human Cytochrome P450 1B1 (CYP1B1) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents fo...
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q16678
Sensitivity	0.063ng/mL
Entrez	1545
Expiration Date	Please enquire.