

Product Datasheet

Human CYP24A1 ELISA Kit (Ready to Use) (orb550445)

Description	Human Cytochrome P450 24A1 (CYP24A1) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents ...
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q07973
Sensitivity	0.053ng/mL
Entrez	1591
Expiration Date	Please enquire.