

Product Datasheet

Human CYP27B1 ELISA Kit (Ready to Use) (orb550448)

Description	Human Cytochrome P450 27B1 (CYP27B1) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents ...
Reactivity	Human
Range	0.78-50ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Sensitivity	0.29ng/mL
Expiration Date	Please enquire.