

Product Datasheet

Human CYP2C18 ELISA Kit (Ready to Use) (orb550450)

Description	Human Cytochrome P450 2C18 (CYP2C18) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents ...
Reactivity	Human
Range	0.625-40ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P33260
Sensitivity	0.247ng/mL
Entrez	1562
Expiration Date	Please enquire.