

Product Datasheet

Human DOCK4 ELISA Kit (Ready to Use) (orb550479)

Description	Human Dedicator Of Cytokinesis 4 (DOCK4) Sandwich ELISA KitReady to Use: Pre-diluted detection reage...
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q8N110
Sensitivity	0.063ng/mL
Entrez	9732
Expiration Date	Please enquire.