

Product Datasheet

Human dSIP ELISA Kit (Ready to Use) (orb550495)

Description	Human Delta-Sleep Inducing Peptide (dSIP) Sandwich ELISA KitReady to Use: Pre-diluted detection reagent...
Reactivity	Human
Range	125-8000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Sensitivity	47pg/mL
Expiration Date	Please enquire.