

Product Datasheet

Human DPYD ELISA Kit (Ready to Use) (orb550534)

Description	Human Dihydropyrimidine Dehydrogenase (DPYD) Sandwich ELISA KitReady to Use: Pre-diluted detection r...
Reactivity	Human
Range	1.56-100ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q12882
Sensitivity	0.65ng/mL
Entrez	1806
Expiration Date	Please enquire.