

Product Datasheet

Human DPAGT1 ELISA Kit (Ready to Use) (orb550557)

Description	Human Dolichyl Phosphate-N-Acetylglucosaminophosphotransferase 1 (DPAGT1) Sandwich ELISA KitReady to...
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q9H3H5
Sensitivity	0.059ng/mL
Entrez	1798
Expiration Date	Please enquire.