

Product Datasheet

Human DUSP9 ELISA Kit (Ready to Use) (orb550573)

Description	Human Dual Specificity Phosphatase 9 (DUSP9) Sandwich ELISA KitReady to Use: Pre-diluted detection r...
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q99956
Sensitivity	0.117ng/mL
Entrez	1852
Expiration Date	Please enquire.