

Product Datasheet

Human EZH1 ELISA Kit (Ready to Use) (orb550633)

Description	Human Enhancer Of Zeste Homolog 1 (EZH1) Sandwich ELISA KitReady to Use: Pre-diluted detection reage...
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q92800
Sensitivity	0.113ng/mL
Entrez	2145
Expiration Date	Please enquire.