

Product Datasheet

Human EPG ELISA Kit (Ready to Use) (orb550652)

Description	Human Epigen (EPG) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for convenient and...
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability tissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q6UW88
Sensitivity	0.112ng/mL
Entrez	255324
Expiration Date	Please enquire.