

Product Datasheet

Human ECRG4 ELISA Kit (Ready to Use) (orb550668)

Description	Human Esophageal Cancer Related Gene 4 (ECRG4) Sandwich ELISA KitReady to Use: Pre-diluted detection...
Reactivity	Human
Range	78.1-5000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q99LS0
Sensitivity	30pg/mL
Entrez	78896
Expiration Date	Please enquire.