

Product Datasheet

Human Estrogen Receptor beta ELISA Kit (Ready to Use) (orb550670)

Description	Human Estrogen Receptor Beta (Estrogen Receptor beta) Sandwich ELISA KitReady to Use: Pre-diluted de...
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q92731
Sensitivity	0.115ng/mL
Entrez	2100
Expiration Date	Please enquire.