

Product Datasheet

Human EAAT2 ELISA Kit (Ready to Use) (orb550677)

Description	Human Excitatory Amino Acid Transporter 2 (EAAT2) Sandwich ELISA KitReady to Use: Pre-diluted detect...
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P43004
Sensitivity	0.113ng/mL
Entrez	6506
Expiration Date	Please enquire.