

Product Datasheet

Human EAAT3 ELISA Kit (Ready to Use) (orb550678)

Description	Human Excitatory Amino Acid Transporter 3 (EAAT3) Sandwich ELISA KitReady to Use: Pre-diluted detect...
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	P43005
Sensitivity	0.117ng/mL
Entrez	6505
Expiration Date	Please enquire.