

## Product Datasheet

### Human FAM20A ELISA Kit (Ready to Use) (orb550694)

|                          |                                                                                                               |
|--------------------------|---------------------------------------------------------------------------------------------------------------|
| <b>Description</b>       | Human Family With Sequence<br>Similarity 20, Member A (FAM20A)<br>Sandwich ELISA KitReady to Use: Pre-dilu... |
| <b>Reactivity</b>        | Human                                                                                                         |
| <b>Range</b>             | 0.156-10ng/mL                                                                                                 |
| <b>Form/Appearance</b>   | Ready to Use                                                                                                  |
| <b>Storage</b>           | All reagents should be stored at -20°C<br>upon receipt                                                        |
| <b>Note</b>              | For research use only                                                                                         |
| <b>Application notes</b> | Pre-diluted detection reagents and<br>enhanced stabilitytissue homogenates<br>or other biological fluids.     |
| <b>Sample Types</b>      | tissue homogenates and other<br>biological fluids.                                                            |
| <b>Uniprot ID</b>        | <b>L8B8N7</b>                                                                                                 |
| <b>Sensitivity</b>       | 0.069ng/mL                                                                                                    |
| <b>Entrez</b>            | <b>54757</b>                                                                                                  |
| <b>Expiration Date</b>   | Please enquire.                                                                                               |