

Product Datasheet

Human FAM3B ELISA Kit (Ready to Use) (orb550695)

Description	Human Family With Sequence Similarity 3, Member B (FAM3B) Sandwich ELISA KitReady to Use: Pre-dilute...
Reactivity	Human
Range	0.625-40ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P58499
Sensitivity	0.235ng/mL
Entrez	54097
Expiration Date	Please enquire.