

Product Datasheet

Human FUCa2 ELISA Kit (Ready to Use) (orb550805)

Description	Human Fucosidase Alpha L2, Plasma (FUCa2) Sandwich ELISA KitReady to Use: Pre-diluted detection reagent...
Reactivity	Human
Range	0.625-40ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q9BTY2
Sensitivity	0.216ng/mL
Entrez	2519
Expiration Date	Please enquire.