

Product Datasheet

Human GPRC5A ELISA Kit (Ready to Use) (orb550819)

Description	Human G Protein Coupled Receptor, Family C, Group 5, Member A (GPRC5A) Sandwich ELISA KitReady to Us...
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q8NFJ5
Sensitivity	0.116ng/mL
Entrez	9052
Expiration Date	Please enquire.