

Product Datasheet

Human GRIA2 ELISA Kit (Ready to Use) (orb550881)

Description	Human Glutamate Receptor, Ionotropic, AMPA 2 (GRIA2) Sandwich ELISA KitReady to Use: Pre-diluted det...
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P42262
Sensitivity	0.115ng/mL
Entrez	2891
Expiration Date	Please enquire.