

Product Datasheet

Human GRIA4 ELISA Kit (Ready to Use) (orb550883)

Description	Human Glutamate Receptor, Ionotropic, AMPA 4 (GRIA4) Sandwich ELISA KitReady to Use: Pre-diluted det...
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	P48058
Sensitivity	0.104ng/mL
Entrez	2893
Expiration Date	Please enquire.