

Product Datasheet

Human GFPT2 ELISA Kit (Ready to Use) (orb550891)

Description	Human Glutamine Fructose-6-Phosphate Transaminase 2 (GFPT2) Sandwich ELISA KitReady to Use: Pre-dilu...
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	O94808
Sensitivity	0.114ng/mL
Entrez	9945
Expiration Date	Please enquire.