

Product Datasheet

Human GPIHBP1 ELISA Kit (Ready to Use) (orb550928)

Description	Human Glycosylphosphatidylinositol Anchored High Density Lipoprotein Binding Protein 1 (GPIHBP1) San...
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q8IV16
Sensitivity	0.054ng/mL
Entrez	338328
Expiration Date	Please enquire.