

Product Datasheet

Human HSPb8 ELISA Kit (Ready to Use) (orb550998)

Description	Human Heat Shock Protein Beta 8 (HSPb8) Sandwich ELISA KitReady to Use: Pre-diluted detection reagen...
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q9UJY1
Sensitivity	0.056ng/mL
Entrez	26353
Expiration Date	Please enquire.