

Product Datasheet

Human HBg2 ELISA Kit (Ready to Use) (orb551010)

Description	Human Hemoglobin Gamma 2 (HBg2) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for c...
Reactivity	Human
Range	312-20000ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, erythrocyte lysates or other biological fluids.
Sample Types	serum, plasma, erythrocyte lysates and other biological fluids.
Uniprot ID	P69892
Sensitivity	95.7ng/mL
Entrez	3048
Expiration Date	Please enquire.