

Product Datasheet

Human HAT1 ELISA Kit (Ready to Use) (orb551054)

Description	Human Histone Acetyltransferase 1 (HAT1) Sandwich ELISA KitReady to Use: Pre-diluted detection reage...
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	014929
Sensitivity	0.103ng/mL
Entrez	8520
Expiration Date	Please enquire.