

Product Datasheet

Human HABP1 ELISA Kit (Ready to Use) (orb551075)

Description	Human Hyaluronan Binding Protein 1 (HABP1) Sandwich ELISA KitReady to Use: Pre-diluted detection rea...
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q07021
Sensitivity	0.072ng/mL
Entrez	708
Expiration Date	Please enquire.