

Product Datasheet

Human ITGb5 ELISA Kit (Ready to Use) (orb551144)

Description	Human Integrin Beta 5 (ITGb5) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for con...
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P18084
Sensitivity	0.065ng/mL
Entrez	3693
Expiration Date	Please enquire.