

## Product Datasheet

### ABCA4 Immunizing peptide (orb55116)

|                          |                                                                            |
|--------------------------|----------------------------------------------------------------------------|
| <b>Description</b>       | ABCA4 Immunizing peptide                                                   |
| <b>Preservatives</b>     | 100µg of Lyophilised peptide, formulation<br>De-ionized (DI) water.        |
| <b>Storage</b>           | Shipped at ambient temperature, store at -<br>20°C                         |
| <b>Note</b>              | For research use only                                                      |
| <b>Application notes</b> | This peptide was used as immunogen to<br>produce orb55116 (ABCA4 antibody) |
| <b>Protein Sequence</b>  | EPECPPQLNTGTQ                                                              |
| <b>Expiration Date</b>   | 12 months from date of receipt.                                            |