

Product Datasheet

Human IL12Rb1 ELISA Kit (Ready to Use) (orb551186)

Description	Human Interleukin 12 Receptor Beta 1 (IL12Rb1) Sandwich ELISA KitReady to Use: Pre-diluted detection...
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P42701
Sensitivity	0.071ng/mL
Entrez	3594
Expiration Date	Please enquire.