

Product Datasheet

Human Leptin Receptor ELISA Kit (Ready to Use) (orb551357)

Description	Human Leptin Receptor (Leptin Receptor) Sandwich ELISA KitReady to Use: Pre-diluted detection reagen...
Reactivity	Human
Range	0.781-50ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P48357
Sensitivity	0.29ng/mL
Entrez	3953
Expiration Date	Please enquire.