

Product Datasheet

Human LNPEP ELISA Kit (Ready to Use) (orb551367)

Description	Human Leucyl/Cystinyl Aminopeptidase (LNPEP) Sandwich ELISA KitReady to Use: Pre-diluted detection r...
Reactivity	Human
Range	1.56-100ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q9UIQ6
Sensitivity	0.58ng/mL
Entrez	4012
Expiration Date	Please enquire.