

Product Datasheet

Human BSDL ELISA Kit (Ready to Use) (orb551382)

Description	Human Lipase, Bile Salt Dependent (BSDL) Sandwich ELISA KitReady to Use: Pre-diluted detection reage...
Reactivity	Human
Range	78.1-5000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	X6R868
Sensitivity	32pg/mL
Entrez	1056
Expiration Date	Please enquire.