

Product Datasheet

Human MGA ELISA Kit (Ready to Use) (orb551458)

Description	Human Maltase Glucoamylase, Intestinal (MGA) Sandwich ELISA KitReady to Use: Pre-diluted detection r...
Reactivity	Human
Range	1.56-100ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q8IWI9
Sensitivity	0.63ng/mL
Entrez	23269
Expiration Date	Please enquire.