

Product Datasheet

Human MPP6 ELISA Kit (Ready to Use) (orb551507)

Description	Human Membrane Protein, Palmitoylated 6 (MPP6) Sandwich ELISA KitReady to Use: Pre-diluted detection...
Reactivity	Human
Range	12.5-800pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q99547
Sensitivity	3.1pg/mL
Entrez	10200
Expiration Date	Please enquire.