

Product Datasheet

Human MYO7B ELISA Kit (Ready to Use) (orb551629)

Description	Human Myosin VIIB (MYO7B) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for conveni...
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q6PIF6
Sensitivity	0.132ng/mL
Entrez	4648
Expiration Date	Please enquire.