

Product Datasheet

Human Nesp1 ELISA Kit (Ready to Use) (orb551676)

Description	Human Nesprin 1 (Nesp1) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for convenien...
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q8NF91
Sensitivity	0.105ng/mL
Entrez	23345
Expiration Date	Please enquire.