

Product Datasheet

Human Neurokinin B ELISA Kit (Ready to Use) (orb551706)

Description	Human Neurokinin B (Neurokinin B) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for...
Reactivity	Human
Range	12.5-800pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q9UHF0
Sensitivity	3.12pg/mL
Entrez	6866
Expiration Date	Please enquire.