

Product Datasheet

Human NKX6-1 ELISA Kit (Ready to Use) (orb551750)

Description	Human NK6 Homeobox Protein 1 (NKX6-1) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents...
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P78426
Sensitivity	0.096ng/mL
Entrez	4825
Expiration Date	Please enquire.