

Product Datasheet

Human NFIB ELISA Kit (Ready to Use) (orb551759)

Description	Human Nuclear Factor I/B (NFIB) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for c...
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	000712
Sensitivity	0.054ng/mL
Entrez	4781
Expiration Date	Please enquire.