

Product Datasheet

Human OSMR ELISA Kit (Ready to Use) (orb551789)

Description	Human Oncostatin M Receptor (OSMR) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents fo...
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityurine, tissue homogenates, cell lysates or other biological fluids.
Sample Types	urine, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q99650
Sensitivity	0.058ng/mL
Entrez	9180
Expiration Date	Please enquire.