

Product Datasheet

Human OSCP1 ELISA Kit (Ready to Use) (orb551793)

Description	Human Organic Solute Carrier Partner 1 (OSCP1) Sandwich ELISA KitReady to Use: Pre-diluted detection...
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q8WVF1
Sensitivity	0.115ng/mL
Entrez	127700
Expiration Date	Please enquire.